

Operational **M**anual

M965/965+ Reader
M965 Mate 2.0 USB
PC Software



Metertech Inc.
Version 2.10

Metertech Inc. provides this publication "as is" without warranty of any kind, either express or implied, including, but not limited to the implied warranties of merchantability or fitness for a particular purpose.

All rights are reserved, including those to reproduce this publication or parts thereof in any form without permission in writing from Metertech Inc.

Information in this publication is subject to change without notice and does not represent a commitment on the part of the vender.

Any errors or omissions which may have occurred in this publication despite the utmost care taken in its production will be corrected as soon as possible, but not necessarily immediately upon detection.

Copyright© Metertech Inc. 2020

Metertech Inc. • 63-2, ChengGong Rd., Sec. 1, NanGang, Taipei, Taiwan, ROC

E-Mail: info@metertech.com • Web Site: www.metertech-inc.com

Tel: 886-2-2783-2854 • Fax: 886-2-2783-1764

Contents

System Requirements	-----	4
Software Installation	-----	5
M965/965+ Instrument Setup	-----	7
To Start M965 Mate 2.0 USB	-----	8
M965 Mate 2.0 USB Menu Software	-----	10
Structure		
Main Window Overview	-----	10
Section A Menu	-----	11
Section B Toolbar	-----	14
Section C Message	-----	14
Section D Temperature monitor	-----	14
Section E Working area & F Data review	-----	15
Section G Special function	-----	15
M965 Mate 2.0 USB Function	-----	16
Basic Parameters	-----	16
Primary and Reference wavelengths	-----	17
Starting method to read plates	-----	17
Plate motion	-----	18
The built-in Incubator	-----	18
The built-in Shaker	-----	18
Well Mapping	-----	19
Quantitative / Cut off / Ratio / Inhibition / QC	-----	20
Print options	-----	24
Main Menu Configuration	-----	25
File menu functions	-----	26
Setup menu functions	-----	29
Help menu functions	-----	33
Toolbar Menu Configuration	-----	35
Message Area Configuration	-----	36
Defining Parameters for Experiment	-----	37
Well Mapping	-----	42
Quantitative Measuring Method	-----	48
Cutoff Measuring Method	-----	51
Ratio/Inhibition Calculation Method	-----	53
Q.C. Calculation Method	-----	54

Printing Options	-----	55
Interpreting the Results	-----	56


System Requirements

- CPU Pentium 4 2.0GHz above for Windows 7 or above.
- 2GB of RAM or above for Windows 7 or above.
- Microsoft .NET Framework 3.5
- 50MB of available hard drive for the program files
- CD ROM drive
- 16bit color display with pixel resolution 1280 x 768 or above.
- Keyboard, Mouse, and RS232 serial port or USB port




Software Installation

To install M965 Mate 2.0 USB


1. Start Windows and close all unnecessary Windows applications.
2. Insert the software CD into the CD-ROM drive. The installer user interface is displayed.
3. Click on the "setup" to install M965 Mate 2.0 USB software.



4. Follow the on-screen instruction.




5. After the installation is completed, click "Close".



M965/965+ Instrument Setup

1. Be sure the M965/965+ instrument is in standalone mode.
Method to switch between standalone and PC modes on M965/M965:
Turn off the instrument first, then press the "OPTION" key while turn on the instrument again, it will switch to other mode.
2. On M965/965+ standalone mode, please go to SETUP / COMPUTER.
3. Move "UP/DOWN" buttons to select USB port and press "ENTER" to confirm selection with mark "S" shown on the right side of that port.
4. Power off and power on the M965/965+ again, then press the "OPTION" key while turn on the instrument, it will switch to PC mode.
5. The instrument starts to do initialization.
6. After initialized, be sure the screen show "pc_usb mode" as figure 1 below.



The image shows a rectangular box representing a screen display. Inside the box, the text is centered and reads: ****ELISA READER **** on the first line, and **pc_usb: off_line** on the second line.

Figure 1

To Start M965 Mate 2.0 USB

Connect the PC and the instrument with an USB cable, then power meter up.

1. Go to M965 mate 2.0 USB


From Start menu →All Programs→Metertech→M965 Mate 2.0 USB




2. For the first time login, key in default value "admin" for both User name and Password in Security window below. Press "OK" to start comport connection.

User name : **admin**

Password : **admin**



3. If the computer appears the message as below, please shut off the M965 mate 2.0 USB software and restart it again, let the software reconnect the reader.



4. If PC software is successfully connected to M965/965+, the "Connected" sign with green background will appear on the upper right message area of the screen.



5. In case the PC software still cannot connect to the instrument, please shut off the PC software and the reader. Be sure the USB cable is plugged in and restart the PC software and Reader again to do the connection.
6. Now, the Reader and software are ready to perform the experiment.

M965 Mate 2.0 USB Menu Software Structure

Main Window Overview

Section A: Menu

Section B: Tool bar


Section C: Message

Section D: Temperature monitor

Section E: Working area

Section F: Data review


Section G: Special function



Section A Menu

File Menu

The File Menu contains the file and print functions for the experiment data and mapping file.



New: Create a new experiment

Load (Data or Map): Load a stored experiment, results or map layout

Save (Data or Map): Save the experiment parameters, results or map layout

Export: Export report to a file with ".csv" file extension


Print: Select experiment to print out

Preview: Preview experiment report format

Exit: Close the M965 Mate 2.0 USB software

Setup Menu

The setup menu contains the M965 Mate 2.0 USB system configuration and user account management.




Security setting: To create a new user ID and set up the security level, or delete the user ID.

Initial: To initial reader and PC software, especially after the reader and the PC software are left for a while, or after re-start the software than to do re-connect.

Help Menu

The help menu provides information on software version, contact information of the



vender, and the user activity records.

About: To provide software version and contact information of the vender.

Audit trail: To record the user activity for trailing.

Section B Toolbar





Plate in/out: To open or close the plate compartment

Shaker: To shake the plate with desired speed and shaking time

Incubator: To control the incubator with desired temperature, and display temperature reading on section D temperature monitor.


Filter tune: To set up the wavelengths of installed filters on 8-slot filter wheel, and have the meter tune the light intensity for each filter.

Post processing: Use the current protocol to re-process data results

Start: To start the experiment with current protocol


Section C Message

During operation, the current status will be shown on the upper right of the screen.




Section D Temperature monitor

When the incubator is activated, the set and actual temperatures are shown on the left of the message area.




Section E Working area & F Data review

The M965 Mate 2.0 USB allows you to define measurement protocols and analyze obtained microplate data. The protocol parameters are input in E Working area, and the test data is shown in F Data review.



Section G Special function

The special function is customized for the Biotest reagent test. This experiment effects only in conjunction with the Biotest reagents.




M965 Mate 2.0 USB Function


Basic Parameters

Measurement types: The M965 Mate 2.0 USB provides three types of measurement, i.e. End point, Two point, and Kinetic measurement.


End point read: During the end point read, the M965/965+ reads at one wavelength, with one-reference wavelength read as optional



Two points read: During the two points read, the M965/965+ reads at two wavelengths, with two-reference wavelength read as optional.



Kinetics read: During the kinetics read, users can define the kinetic method by selecting Average rate, Maximum rate, Maximum Abs, Total delta Abs, Time to max rate, or Time to max Abs in the Kinetic method list. The user can also define the measure numbers and interval.



To Set up a measurement with End point, Two points or Kinetic method, user need to define following parameters.

Primary and Reference wavelengths

If a Primary wavelength is defined alone, the M965/965+ reads the plate only once at a single wavelength. If a Reference wavelength is defined, the plate will be read twice and automatically calculate the delta Abs between these two readings.

Method to set up the Primary and Reference wavelengths:

1. Select the Measurement type of End Point, or Two Points.
2. Enter the Primary wavelength in Main1 or Main2, and the reference Ref1 or Ref2

Starting method to read plates

If the "Immediate" option is selected, the instrument starts reading the plate right after the Start button is pressed. Users can also define the period of the plate reading delay.

To define the starting method,

1. Select the "Immediate" option.
2. Or choose "Delay", then input the delay time in second.

Plate motion

Users can select the plate motion as stepping in milliseconds or continuous mode.

The built-in Incubator

The incubator will keep the plate stay at temperature-controlled environment.

Users can activate the incubator by

1. Clicking the incubator button on toolbar to display the incubator pop-up menu.
2. Entering the desired temperature on the pop-up menu, and press "Activate" tab to start the temperature control.

The built-in Shaker

The built-in Shaker in the instrument allow user to define speed setting as Low 8Hz, Medium 11Hz, or High 14Hz. Users can also define the shaking period.


To enable the shaker,

1. Click the shaker button
2. Select the speed to be Low, Medium or High
3. Define the shaking period in second.

Well Mapping

Users can define five types of different wells. They are Blank, Standard, Sample, Positive, and Negative on the Type menu at Map layout.

Users can also define the sample well name at Map layout as required.



Map layout

	1	2	3	4	5	6	7	8	9	10	11	12
A	SAM01	SAM02	SAM03	SAM04	SAM05	SAM06	SAM07	SAM08	SAM09	SAM10	SAM11	SAM12
B	SAM01	T01-2 *										
C	SAM01	T01-3 *										
D				POS01	STD01	STD02	STD03	STD04	STD05			

test-A

Map layout


	1	2	3	4	5	6	7	8	9	10	11	12
A	SAM01	SAM02	SAM03	SAM04	SAM05	SAM06	SAM07	SAM08	SAM09	SAM10	SAM11	SAM12
B	SAM01	T01-1 *										
C	SAM01	T02-3										
D				POS01	STD01	STD02	STD03	STD04	STD05			

Quantitative / Cut off / Ratio / Inhibition / QC

Quantitative setting

The M965 Mate 2.0 USB allows user to define quantitative analysis to determine the sample concentration. Seven types of curve fitting equations are built to calculate standard polynomial coefficients. Users can select Curve on plate, Stored curve, Standard line, or Concentration factor to define Quantitative method.


1. Click the Quantitative / Cut off / Ratio / Inhibition / QC button.
2. Click the check mark in front of the Quantitative settings area, and be sure the check mark turned into green.
3. Define the desired parameters.



Cutoffs

Cutoffs are used to classify results. Users can define three Cutoff methods as Single, Double, Calculation and S/CO ratio


1. Click the Quantitative / Cut off / Ratio / Inhibition / QC button.
2. Click the check mark in front of the Cutoff settings area, and be sure the check mark turned to green.
3. Define the desired parameters.



Ratio/Inhibition

The M965 Mate 2.0 USB will take a reference (B0) and other samples (B) to calculate the Ratio/Inhibition


1. Click the Quantitative / Cut off / Ratio / Inhibition / QC button.
2. Click the check mark in form of the Ratio/Inhibition settings area, and be sure the check mark turned green.
3. Define the desired parameters.



Q.C.


The M965 Mate 2.0 USB provides Q.C. algorithm for experiment to determine the results.

1. Click the Quantitative / Cut off / Ratio / Inhibition / QC button.
2. Click the check mark in front of the Q.C. settings area, and be sure the check mark turned green.
3. Define the desired parameters.






Print options

Users can define the Project title, User name, experiment Note, and check the desired items in Sections for printing to print the result of the experiment.



Main Menu Configuration


There are three functions on the main menu. They are File, Setup and Help.




File menu functions

There are seven options, i.e. New, Load (Data or Map), Save (Data or Map), Export, Print, Preview, and Exit under the main menu.


1. New: Create a new experiment
2. Load (Data or Map) : Load an existing experiment file or map layout



3. Save (Data or Map) : Save experiment file or map layout



4. Export: To export to a file with file extension "csv", and it can be loaded into the excel, notepad or google spreadsheet
5. Print: To print report using the printer connected to the PC
6. Preview: To preview the experiment report before printing



7. Exit: to end the M965 Mate 2.0 USB operation



Setup menu functions

The setup menu includes the Security setting and the Initial.


1. Security setting: Users should log in with their own ID and Password to start the experiment. This function lets you log in / log out system, create, and or delete user ID.

Logo in system

To enter the system, please use the default value (admin) to log in.

User name: admin

Password: admin




A screenshot of a Windows-style dialog box titled "Security". It contains a "Login" section with three input fields: "User name" containing "admin", "Password" containing "*****", and "Project name" which is empty. An "OK" button is located at the bottom right of the dialog.

Entering the wrong user name or password, The M965 Mate 2.0 USB will show the pop-up window as below.




Log out system

To log out system or change the operator ID




Create a new ID

This function can be performed by users with administrator authority only.



Enter the new ID with user name, password, and ID authority (administrator or user), and then press Save button to accept.




The image shows a 'Security' dialog box with a title bar containing a lock icon and a close button. The main area is titled 'New user' and contains three text input fields: 'User name : Test 01', 'Password : *****', and 'Confirm password : *****'. Below these fields are two radio buttons: 'Administrator' (selected) and 'User'. At the bottom right of the dialog are two buttons: 'Cancel' and 'Save'.

Successfully creating a new user will see the pop-up window as below.




Delete ID

This function can be performed by users with administrator authority only.




Enter the ID which you want to delete and press Delete to execute.



2. Initial

This function is re-connect the M965/965+ Reader and the computer.

Especially after the reader and the PC software are left for a while, or after re-start the software than to do re-connect.




Help menu functions

The help menu, not only shows the information about the vendor and PC software version, but also trails the log-in and operation record.

1. About: Provide the contact information of vendor and PC software version




2. Audit trail: To review the records of activities of different ID




Toolbar Menu Configuration

There are Plate In/Out, Shaker, Incubator, Filter tune, and Post processing tabs on the toolbar menu.


1. Plate in/out: To move the plate holder in or out, the plate holder status will show on the status bar



2. Shaker: This tab is used to configure and operate the shaker. The shaker has three translation speeds, i.e. low (8Hz), Medium(11Hz) and High (14Hz)




3. Incubator: To warm up the incubator at set-point temperature from the lowest 15°C to the highest 50°C. If the ambient temperature is higher than 15°C, the effective lowest temperature should be set to the ambient temperature + 3C.



4. Filter tune: The 965Mate has an eight- slot filter wheel for user to install filters



5. Post Processing: Use the current parameters of selected measure mode and recalculate the data



Message Area Configuration


The message area contains two parts, the status message and the temperature monitor.

1. Status message: To display the status of instrument current operation condition. All messages are listed in the following chart.



Message	Description
Initializing	Initializes the instrument
Standby	Test ready
Data reading	Load data
Post processing	Recalculate the data
Disconnect	The instrument has no connection with M965 Mate 2.0 USB
Connected	The instrument connects with M965 Mate 2.0 USB
Filter tune	Start tuning filter

2. Temperature monitor: To display the real-time temperature and set-point temperature of instrument incubator.




Defining Parameters for Experiment


Defining Parameters

When starting an experiment, users must first define the parameters such as wavelength, reading method, plate motion, incubator, and shaking. Above functions are included in Basic parameters tab.


1. Measurement Type: Users can define three measuring types, i.e. End point, Two point and Kinetic.




a、 End Point




b、 Two Points




c、 Kinetics: Kinetics measuring method can only select main filter without reference filter




2. Filter wavelength: Users need to select the main filter wavelength for the desired experiment. In addition, users can also select a reference wavelength.



3. Starting method: Define when to start the selected experiment.
 - a、 Immediate: Start measurement right after pressing the START tab
 - b、 Delay: Users can define 0~999s as delay time before starting measurement.




4. Plate motion: To define how the plate is moved when measuring
- a. Continuous: When measuring, the plate is translated smoothly during the entire motion stroke.





- b. Stepping: User can define the stepping intervals among 0~999 msec. In kinetic mode, there are variable stepping interval and fixed stepping interval to be selected.




5. Incubation: Users can define the incubator temperature by clicking the incubation tab. The temperature can be set from ambient 15°C to 50°C.




6. Shaker: The shaker of the instrument shakes with three types of speed, and the shaking time can be arranged among 0~999s



7. Two point interval: Users can select the two point interval among 3~999s



8. Kinetic method, numbers, and interval: In kinetic measurement mode, user can select the data calculation method, test cycles, and cycle interval.
- a. Kinetic method: Users can select Average rate, Maximum rate, Maximum OD, Total delta OD, Time to max slope, Time to max OD for mapped wells calculation.



- b、 Measure number: User can enter the measuring numbers among 3~255 cycles.
- c、 Measure interval: User can enter the measure interval. They are among 4~65535s in continuous motion, and 6~65535s in stepping motion.


Kinetic method	Average rate
Measure numbers	3
Measure interval	4 sec.

Well Mapping

The M965 Mate 2.0 USB provides five types of well for the user to define 96-well map, also including define each sample well ID. Moreover, the user can save the mapped wells and reload them for further uses.

1. Save and load map layout:


Users can load or save their map layout from File/ Load (Data or Map) or Save (Data or Map) functions



2. Well mapping method:

- a. Select the well type to be defined (Blank, Positive, Negative, Sample, Standard) on the map layout. On the right side, enter the concentration values if standard is selected.

Note: Sample is the only available type in kinetic mode.



- b. Select the fill and replicate directions, enter next fill number, and replicate number.

Fill direction: Row Column

Replicate direction: Row Column

- c. Use mouse to draw an area, which wells are to be placed with selected type.
- d. Right click on the mouse to select the fill option.

Map layout

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

Type: Standard
 Fill direction: Row Column
 Replicate direction: Row Column
 Next fill number: 1
 Replicate number: 1

STD Conc.
 STD01 N.A.
 STD02 N.A.
 STD03 N.A.
 STD04 N.A.
 STD05 N.A.
 STD06 N.A.
 STD07 N.A.
 STD08 N.A.
 STD09 N.A.
 STD10 N.A.
 STD11 N.A.
 STD12 N.A.
 STD13 N.A.
 STD14 N.A.
 STD15 N.A.
 Input range: 0.001-99999.9

- e. The selected ten standards are thus located on the well map.

Map layout



	1	2	3	4	5	6	7	8	9	10	11	12
A	STD01 C01-1	STD06 C06-1										
B	STD02 C02-1	STD07 C07-1										
C	STD03 C03-1	STD08 C08-1										
D	STD04 C04-1	STD09 C09-1										
E	STD05 C05-1	STD10 C10-1										
F												
G												
H												

Type: Standard
 Fill direction: Row Column
 Replicate direction: Row Column
 Next fill number: 11
 Replicate number: 1

STD Conc.
 STD01 1.000
 STD02 1.000
 STD03 1.000
 STD04 1.000
 STD05 1.000
 STD06 1.000
 STD07 1.000
 STD08 1.000
 STD09 1.000
 STD10 1.000
 STD11 N.A.
 STD12 N.A.
 STD13 N.A.
 STD14 N.A.
 STD15 N.A.
 Input range: 0.001-99999.9

3. Fill and replicate well:
 - a. Fill direction: To number the sequence of selected sample type in column or row direction.
 - b. Replicate direction: To number the sequence of the replicates of selected sample type in column or row direction.

Example of filling and replicating the well map



4. Blank, positive control, and negative control each has only one name (BLK01, POS01, NEG01).
5. Standard can be configured as 1~15 names(STD01~STD15)
6. Sample has 96 defaults names most (SAM01~SAM96).

7. Types of well:

- a、 BLK: Which is painted with light green background on the well map
 - b、 POS: Which is painted with light red background on the well map
 - c、 NEG: Which is painted with light blue background on the well map
 - d、 Sample: Which is painted with light orange background on the well map
 - e、 Standard: Which is painted with light purple background on the well map
- Users must fill in the concentration values of selected standards in ascending or descending order.

Map layout

	1	2	3	4	5	6	7	8	9	10	11	12
A	BLK01 Z01-1	BLK01 Z01-2	BLK01 Z01-3	BLK01 Z01-4	BLK01 Z01-5							
B	POS01 P01-1	POS01 P01-2	POS01 P01-3	POS01 P01-4	POS01 P01-5							
C	NEG01 N01-1	NEG01 N01-2	NEG01 N01-3	NEG01 N01-4	NEG01 N01-5							
D	SAM01 T01-1	SAM02 T02-1	SAM03 T03-1	SAM04 T04-1	SAM05 T05-1							
E	STD01 C01-1	STD02 C02-1	STD03 C03-1	STD04 C04-1	STD05 C05-1							
F												
G												
H												

STD Conc.

STD01	1.000
STD02	1.000
STD03	1.000
STD04	1.000
STD05	1.000
STD06	N.A
STD07	N.A
STD08	N.A
STD09	N.A
STD10	N.A
STD11	N.A
STD12	N.A
STD13	N.A
STD14	N.A
STD15	N.A

Type: Fill direction: Row Column Replicate direction: Row Column Next fill number: Replicate number: Input range: 0.001-99999.9


8. Edit sample ID

- a. After creating Sample wells at the Map layout, right click the mouse to select the Edit ID option at the pop-up window, and key in the sample ID as required.

	1	2	3	4	5	6	7	8	9	10	11	12
A	SAM01 T01-1	SAM07 T07-1	SAM01 T01-2	SAM02 T02-2	SAM03 T03-2							
B	SAM02 T02-1	SAM08 T08-1	SAM01 T01-1	SAM02 T02-3								
C	SAM03 T03-1											
D	SAM04 T04-1											
E	SAM05 T05-1											
F	SAM06 T06-1	SAM07 T07-2	SAM08 T08-2	SAM09 T09-1	SAM10 T10-1							
G												
H												

Context menu options: Fill, Modify, Clear, Clear Group, Clear All, Edit ID

	1	2	3	4	5	6	7	8	9	10	11	12
A	SAM01 T01-1	SAM07	SAM01	SAM02	SAM03							
B	SAM02 T02-1											
C	SAM03 T03-1											
D	SAM04 T04-1	T10-1										
E	SAM05 T05-1	SAM11 T11-1										
F	SAM06 T06-1	SAM12 T12-1										
G												
H												



- b. Check the sample ID: Those named sample wells will come up a * mark behind the replicate number. Moving the mouse cursor to the well and the ID will pop up automatically.

	1	2	3	4	5	6	7	8	9	10	11	12
A	SAM01 T01-1 *	SAM07 T07-1	SAM01 T01-2 *	SAM02 T02-2	SAM03 T03-2							
B	SAM01 T02-1 *	SAM01 T01-1 *	SAM01 T01-3 *	SAM02 T02-3								
C	SAM01 T03-1	T09-1	SAM01 T01-4 *									
D	SAM04 T04-1	SAM10 T10-1										
E	SAM05 T05-1	SAM11 T11-1										
F	SAM06 T06-1	SAM12 T12-1										
G												
H												

- c. Auto ID naming of replicated wells: If the sample has several replicates, it is only to name any picked one, and the rest ID will be automatically created.

	1	2	3	4	5	6	7	8	9	10	11	12
A	SAM01 T01-1 *	SAM07 T07-1	SAM01 T01-2 *	SAM02 T02-2	SAM03 T03-2							
B	SAM02 T02-1	SAM08 T08-1	SAM01 T01-3 *	SAM01 T01-2 *								
C	SAM03 T03-1	SAM09 T09-1	SAM01 T01-4 *	[2018061110]								
D	SAM04 T04-1	SAM10 T10-1										
E	SAM05 T05-1	SAM11 T11-1										
F	SAM06 T06-1	SAM12 T12-1										
G												
H												

•	1	2	3	4	5	6	7	8	9	10	11	12
A	SAM01 T01-1 *	SAM07 T07-1	SAM01 T01-2 *	SAM02 T02-2	SAM03 T03-2							
B	SAM02 T02-1	SAM08 T08-1	SAM01 T01-3 *	SAM02 T02-3								
C	SAM03 T03-1	SAM09 T09-1	SAM01 T01-4 *									
D	SAM04 T04-1	SAM10 T10-1										
E	SAM05 T05-1	SAM11 T11-1										
F	SAM06 T06-1	SAM12 T12-1										
G												
H												

- d. Amend or cancel sample ID: Right click the mouse to select the Edit ID option, then amend the ID or remove ID.

•	1	2	3	4	5	6	7	8	9	10	11	12
A	SAM01 T01-1 *	SAM07 T07-1	SAM01 T01-2 *	SAM02 T02-2	SAM03 T03-2							
B	SAM02 T02-1	SAM08 T08-1	SAM01 T01-3 *									
C	SAM03 T03-1	SAM09 T09-1	SAM01 T01-4 *									
D	SAM04 T04-1	SAM10 T10-1										
E	SAM05 T05-1	SAM11 T11-1										
F	SAM06 T06-1	SAM12 T12-1										
G												
H												

ID name

Test

OK


•	1	2	3	4	5	6	7	8	9	10	11	12
A	SAM01 T01-1 *	SAM07 T07-1	SAM01 T01-2 *	SAM02 T02-2	SAM03 T03-2							
B	SAM02 T02-1	SAM08 T08-1	SAM01 T01-3 *	SAM02 T02-3								
C	SAM03 T03-1	SAM09 T09-1	SAM01 T01-4 *									
D	SAM04 T04-1	SAM10 T10-1										
E	SAM05 T05-1	SAM11 T11-1										
F	SAM06 T06-1	SAM12 T12-1										
G												
H												

[Test]


Quantitative Measuring Method

The M965 Mate 2.0 USB provides four types of Quantitative method, i.e. Curve on plate, Stored curve, Standard line, and Concentration factor.

1. Curve on plate: Use the standard on the well plate for the calibration curve calculation. There are seven types of curve fitting equations on the M965 Mate 2.0 USB
 - a、 Linear regression
 - b、 Quadratic polynomial
 - c、 Cubic polynomial
 - d、 Point to point
 - e、 Cubic spline
 - f、 2 parameters logit-log
 - g、 4 parameters logistic
 - h、 Advance BIO (Option function)




2. Stored curve: Users can load their stored curve for quantitative measurement, these curves with file extension ".cuv".



3. Standard line: User can use the $Abs=A* Conc+B$ equation, and enter the values of A and B to calculate a standard line.


The value of A can be : -999999.999 ~ +999999.999

The value of B can be : -999999.999 ~ +999999.999



4. Concentration factor: User can enter a factor for calculating the concentration.

The value of F can be : -999999.999 ~ +999999.999



5. Measurement unit: Users can select 15 types of measurement unit "None" , "G/dL" , "U/L" , "G/L" , "ug/dL" , "ABS" , "mg/dL" , "OD" , "mABS" , "U/mL" , "ug/mL" , "mEq/L" , "mmol/L" , "umol/L" , "ng/mL". When "None" is selected, user can enter the desired measurement unit

concentration unit

Unit

User key in


Cutoff Measuring Method

The M965 Mate 2.0 USB provides four types of Cutoff measuring method.

1. Single cutoff method: User can enter a threshold of 0.0000~4.0000, and define OD result to be positive or negative.



2. Double cutoff method: Users can define the high and low thresholds. The high and low values can be among 0.0000~4.0000. The M965 Mate 2.0 USB determines OD results that are higher than high threshold, lower than low threshold, or between low and high thresholds to be positive(+), negative(-), or in-between (*) respectively.



3. Calculation cutoff method: User can create a maximum of four formulas as the thresholds calculation and categorize the OD readings into 5 groups.

The equation listed below is applied to construct the thresholds with given a, b and c values.

$$EQ_n = a * PC + b * NC + c,$$


Where PC means Positive Control, and NC means Negative Control.

The value of a, b and c can be -1000.000 ~ +1000.000

The calculated threshold values must follow the rule below:

$$EQ1 > EQ2 > EQ3 > EQ4$$

Example: With four thresholds applied, the OD reading higher than EQ1, between EQ1 and EQ2, between EQ2 and EQ3, between EQ3 and EQ4, or below EQ4 is labeled by "++", "+", "* ", "-" or "--" respectively.




4. S/CO ratio method: User uses the ratio of sample OD divided by CO value as the benchmark, and also key in the high and low thresholds. The M965 Mate 2.0 USB will then determine S/CO ratio that are higher than high threshold, lower than low threshold, or between low and high thresholds to be positive(+), negative(-), or in-between (*) respectively.

The equation listed below is applied to construct the CO value with given a, b and c values.

$$\boxed{CO} = a * \boxed{PC} + b * \boxed{NC} + c$$

Where \boxed{PC} means Positive Control, and \boxed{NC} means Negative Control.

The value of a, b and c can be -1000.000 ~ +1000.000




Ratio/Inhibition Calculation Method

Select B0 as the standard value to calculate the rest of the plate well Bn

1. Ratio/Inhibition operating procedure

- a、 Ratio = $(B_n/B_0)\%$
- b、 Inhibition = $100\% - (B_n/B_0)\%$
- c、 Must have sample on B0 position or the software will show error
- d、 If the selected B0 has replicate number greater than one, the actual B0 value will be the average reading of this sample.
- e、 If B0 value is 0, the software will show error
- f、 If ratio is over 200%, the software will show HI; if lower than -200%, the software will show LO



Q.C. Calculation Method

The purpose of the QC Calculation Method is to determine the reliability of the experiment.

1. At most 4 equations are applied to obtain the calculation results, QC1, QC2, QC3 and QC4.
2. Combining above QCs with logic operators OR, AND, and XOR to obtain the QC calculation result. The truth or falseness of QC decides the experiment to be pass or fail.
3. The value of a can be -1000.000 ~ +1000.000
4. The value of b can be -1000.000 ~ +1000.000
5. The value of c can be 1000.000 ~ +1000.000
6. The value of H can be -9999999.999 ~ +9999999.999
7. The value of L can be -9999999.999 ~ +9999999.999

The screenshot shows the software interface for the Q.C. Calculation Method. The interface is organized into several sections:

- Protocol select function:** Includes buttons for 'Basic parameters', 'Quantitative / Cut off / Ratio/Inhibition / Q.C.' (highlighted with a red box), 'Print options', 'Data analysis', and 'Special function'.
- Quantitative settings:** Includes 'Curve on plate' (highlighted with a green checkmark), 'Stored curve', 'Standard line', and 'Concentration factor'. It also has options for 'Data curve fit' (Linear regression), 'X-scale' (Linear, Log10), and 'Y-scale' (Linear, Log10).
- Cutoff settings:** Includes 'Single', 'Double', 'Calculation', and '(S/CO) ratio' options. It also has a 'Threshold' field and options for 'Cutoff label' (+/-).
- Ratio / Inhibition settings:** Includes 'Ratio: B/B0 %' and 'Inhibition: 100-B/B0 %' options, along with a 'Definition of B0' field.
- Q.C. settings:** Includes a 'General equation' field with the formula $L \leftarrow aPC + bNC + c \leftarrow H$. It also has a table for setting up QCs (QC1, QC2, QC3, QC4) with checkboxes and dropdown menus for each parameter (L, a, b, c, H).

The 'START' button is highlighted with a green box. The status bar at the top right shows 'Standby' and 'Connected'.

Printing Options

Users can input project name, operator name, and experiment note to differentiate experiment reports. Users can also check boxes in the Section for printing to determine which items need be printed on the report.

File Setup Help

Plate in/out Shaker Incubator Filter tune Post processing

Protocol select function

Basic parameters

Quantitative / Cut off / Ratio/Inhibition / Q.C.

Print options

Data analysis

Data

Special function

Advance BIO

START

Print options

Title setting
date/time of the performed measurement :

Project name :
User : admin
Note :

Sections for printing


Title Data Q.C.
 Results Cut off Kinetics
 Calibration Ratio/Inhibition

Select all
Cancel all


Interpreting the Results

The M965 Mate 2.0 USB will generate the result data after the experiment is completed. Press the Data tab on the left window, and select tab Results, Calibration, Data, Cutoff, Ratio/ Inhibition, Q.C, or Kinetic to view their experiment results.

1. Results: Click on the Results tab to review the parameter setup, plate layout, Raw OD, and Con Matrix of the experiment.




2. Calibration: When Quantitative is checked, calibration curve will be displayed according to the setting parameters.



- a. Layout: Shows the well mapping layout of the plate. Different types of well uses a different color to represent.
- b. Source data: Shows the source data for the quantitative measurement.

- i. In end point measurement, if there is no reference filter then the main filter (M1) data is the source data. If there is reference filter then M1 – R1 is the source data.
 - ii. In Two points measurement if there is no reference filter, the source data will be M1
 - iii. In Two points measurement if there is reference filter then the source data will be D1=M1-R1
 - iv. During Kinetic measurement, user cannot use reference filter, the M1 data will be the source data
- c. Calibrators: Use C01~C15 to represent each STD's name and OD value, and show the average measurement and the standard concentration value.
 - d. Calib Curve: When using standard curve (Curve on plate or stored curve), apply selected fitting method to create a standard curve and its coefficients.
 - e. Residuals table: Use C01~C15 to show standard concentration values (C set), Average Abs, calculated concentration (Ccal), and their difference (Ccal-Cset).
 - f. Curve Viewer: User can double click on the curve to enable the curve viewer. User can also store the curve by pressing the Save Curve tab on the right. The default curve is stored in M965 Mate 2.0 USB \StdCurve directory.



3. Data sheet: The raw data and calculated results of entire mapped wells can be listed in one data sheet. The sheet provides information about Name, Well location, Replicate numbers, Abs, SD, CV%, Conc, Measuring unit, Cutoff, Inhibition % and Well ID name. The average of replicated data is displayed by "_avg" next to the well ID.

Name	well	Replicate	Abs.	SD	CV%	Conc.	Unit	Cutoff	RI(%)	ID
POS. CONT.										
POS1	D4	1	0.009	---	---	---	---	---	---	---
POS1	E4	2	0.001	---	---	---	---	---	---	---
POS1_avg	---	---	0.005	0.004	80.00	212.964	None	N/A	N/A	---
NEG. CONT.										
NEG1	G4	1	0.002	---	---	---	---	---	---	---
NEG1	H4	2	0.018	---	---	---	---	---	---	---
NEG1_avg	---	---	0.010	0.006	80.00	215.946	None	N/A	N/A	---
SAMPLES										
SAM1	A1	1	-0.217	---	---	---	---	---	---	---
SAM1	B1	2	-0.217	---	---	---	---	---	---	---
SAM1	C1	3	-0.217	---	---	---	---	---	---	---
SAM1_avg	---	---	-0.217	0.000	ERR	-57.726	None	N/A	N/A	---
SAM2	A2	1	-0.217	---	---	---	---	---	---	---
SAM2	B2	2	-0.217	---	---	---	---	---	---	---
SAM2	C2	3	-0.217	---	---	---	---	---	---	---
SAM2_avg	---	---	-0.217	0.000	ERR	-57.726	None	N/A	N/A	---
SAM3	A3	1	-0.217	---	---	---	---	---	---	---
SAM3	B3	2	-0.217	---	---	---	---	---	---	---
SAM3	C3	3	-0.217	---	---	---	---	---	---	---
SAM3_avg	---	---	-0.217	0.000	ERR	-57.726	None	N/A	N/A	---
SAM4	A4	1	0.015	---	---	---	---	---	---	---

4. Cutoff results: Clicking the Cutoff tab, the M965 Mate 2.0 USB shows the cutoff symbols on mapped wells. Depending on the conditions, there will be five symbols to represent the cutoff results.


{ ++ } 、 { + } 、 { * } 、 { - } 、 { -- }

	1	2	3	4	5	6	7	8	9	10	11	12
A	-	-	-	-	-	-	+	+	+	-	-	-
B	-	-	-	-	-	-	+	+	+	-	-	-
C	-	-	-	-	-	-	+	+	+	-	-	-
D							+	+	+			
E							-	+	+			
F							-	+	+			
G							-	+	+			
H							-	+	+			


5. Ratio/Inhibition results: Clicking the Ratio/Inhibition tab, the M965 Mate 2.0 USB shows ratio or inhibition values of mapped wells. Data higher than 200% is shown Hi, and lower than -200% is shown LO.

	1	2	3	4	5	6	7	8	9	10	11	12
A	NaN%	NaN%	NaN%	-Infinity%	-Infinity%	-Infinity%	-Infinity%	-Infinity%	-Infinity%	NaN%	NaN%	-Infinity%
B	NaN%	NaN%	NaN%	-Infinity%	-Infinity%	-Infinity%	-Infinity%	-Infinity%	-Infinity%	NaN%	NaN%	NaN%
C	NaN%	NaN%	NaN%	-Infinity%	-Infinity%	-Infinity%	-Infinity%	-Infinity%	-Infinity%	NaN%	NaN%	NaN%
D				-Infinity%	-Infinity%	-Infinity%	-Infinity%	-Infinity%	-Infinity%			
E				-Infinity%	-Infinity%	-Infinity%	-Infinity%	-Infinity%	-Infinity%			
F				-Infinity%	-Infinity%	-Infinity%	-Infinity%	-Infinity%	-Infinity%			
G				-Infinity%	-Infinity%	-Infinity%	-Infinity%	-Infinity%	-Infinity%			
H				-Infinity%	-Infinity%	-Infinity%	-Infinity%	-Infinity%	-Infinity%			

- Q.C results: Clicking the QC calculation method, the M965 Mate 2.0 USB shows the QC criteria, Pass condition, and Result on the data sheet.



- Kinetic results: When using the kinetic measuring method, M965 Mate 2.0 USB will display the kinetic curves for each mapped wells. User can check the reaction rate easily on this screen.



Double click on the selected well to show a detailed view of the well number and OD value at selected sampling number.

